



RHEA MEDIUM

\$2,349.00

The Rhea features an open-arm design that holds numerous features/accessories, MaxiComfort® Positioning and HeatWave™ Technology. Available in two sizes.

SKU: PR442-MED

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Zero Gravity \(MaxiComfort\)](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medical Device	Yes
Weight Capacity	375
Positioning	MaxiComfort
Comfort Zone(s)	4
Twilight Tilt Technology	NO
ZG+	NO
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	NO
Arm Rest Covers	No



Overall Width	37"
Overall Height	45"
Seat Width	22"
Floor to Top of Seat	20"
Seat Depth	21"
Seat to Top of Back	29"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Right Arm Compartment
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	NO
HeatWave Technology Option	Standard
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	NO
Straight Lift Option	Standard



USB Charging Port	Left Arm Compartment
Battery Backup	Lithium-Ion
SmartTek Diagnostic System	Standard
Place of Origin	Imported