

FDA Class II Medical Device



RHEA WITH TWILIGHT & NIRVANA

\$3,259.00 - \$3,379.00 Price range:
\$3,259.00 through \$3,379.00

Total-body wellness made effortless. The Rhea features patented Twilight Tilt positioning, Nirvana™ Swedish Massage with Far Infrared Heat, plus built-in accessories.

SKU: PR447

Categories: [5 Zone Chairs](#), [DeLuna](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

VARIATIONS

FDA Class II Medical Device



PR447-LAR \$3,379.00

Available in 2 sizes, the Rhea features our patented Twilight Tilt Technology & MaxiComfort® Positioning, Nirvana Swedish Massage & Far Infrared HeatWave Technology.

Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
	PR447-MED	\$3,259.00			Available in 2 sizes, the Rhea features our patented Twilight Tilt Technology & MaxiComfort® Positioning, Nirvana Swedish Massage & Far Infrared HeatWave Technology.	Medium

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Large, Medium
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	Twilight Tilt
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	37"
Overall Height	45"

Seat Width	22"
Floor to Top of Seat	20"
Seat Depth	21"
Seat to Top of Back	29"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Right Arm Compartment
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	No
Nirvana with HeatWave Option	Standard
HeatWave Technology Option	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard

USB Charging Port	YES
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	Imported