

FDA Class II Medical Device



RHEA

\$2,449.00 - \$2,549.00 Price range:
\$2,449.00 through \$2,549.00

The Rhea features an open-arm design that holds numerous features/accessories, MaxiComfort® Positioning and HeatWave™ Technology. Available in two sizes.

SKU: PR442

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [New](#)

VARIATIONS

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<p>FDA Class II Medical Device</p> 	PR442-LAR	\$2,549.00			<p>The Rhea features an open-arm design that holds numerous features/accessories, MaxiComfort® Positioning and HeatWave™ Technology. Available in two sizes.</p>	Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<small>FDA Class II Medical Device</small> 	PR442-MED	\$2,449.00			The Rhea features an open-arm design that holds numerous features/accessories, MaxiComfort® Positioning and HeatWave™ Technology. Available in two sizes.	Medium

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Large, Medium
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	4
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	37"
Overall Height	45"



Seat Width	22"
Floor to Top of Seat	20"
Seat Depth	21"
Seat to Top of Back	29"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Right Arm Compartment
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	No
HeatWave Technology Option	Standard
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	Left Arm Compartment
Battery Backup	Lithium-Ion

**SmartTek Diagnostic
System**

Standard

Place of Origin

Imported