

FDA Class II Medical Device



TITAN WITH TWILIGHT MEDIUM

\$3,299.00

The Titan offers our patented Twilight Tilt Technology, MaxiComfort® Positioning and additional accessories for an enhanced relaxation experience.

SKU: PR449-MED

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medical Device	Yes
Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No

Overall Width	38.5"
Overall Height	46"
Seat Width	20.5"
Floor to Top of Seat	20.5"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	No
HeatWave Technology Option	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard

USB Charging Port	Outside Arm
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	Imported