



TITAN MEDIUM

\$2,349.00

The Titan offers our MaxiComfort® Positioning and numerous accessories to provide hours of unmatched comfort and convenience.

SKU: PR448

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [Zero Gravity](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medial Device	YES
Weight Capacity	375
Positioning	Twilight
Comfort Zone(s)	4
Twilight Tilt Technology	NO
ZG+	NO
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	NO
Arm Rest Covers	No



Overall Width	38.5"
Overall Height	46"
Seat Width	20.5"
Floor to Top of Seat	20.5"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Arm
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	NO
HeatWave Technology Option	NO
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	NO
Straight Lift Option	Standard



USB Charging Port	Outside Arm
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	Imported