

FDA Class II Medical Device



AVENGER

\$4,699.00

The Golden Avenger—a four-wheel, heavy-duty scooter—provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature's most challenging terrains.

SKU: GA541

Categories: [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range	18 miles
Ground Clearance	4"
Turning Radius	67"
Type Batteries	NF22
Colors	Red, Silver
Features	Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, PGS 140 Amp Controller
Charger	Off Board



Length	57.5"
Width	24.5"
Height (ground to top of back of seat)	42.5" - 45.5"
Height (ground to top of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
Number Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	38.5"
Weight of Unit (assembled)	265 lbs
Front Half Weight	68 lbs
Rear Half Weight	68 lbs
Seat with Arms Weight	55 lbs
Batteries Weight	37 lbs ea
Heaviest Piece	68 lbs
Tires (standard)	Pneumatic
Front Tires	13"
Rear Tires	13"
Casters	N/A
Rear Anti-Tippers	3"
Standard Seat	Captain's Seat
Seat Back Height - no headrest	21" - 24"
Back Height with Headrest	27.5" - 31.5"
Seat Width x Depth	20" X 17"
Seat Color	Gray Vinyl
Warranty Drive Train	2 years



Warranty Batteries

1 year