



## AVENGER

**\$4,699.00**

The Golden Avenger—a four-wheel, heavy-duty scooter—provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature’s most challenging terrains.

---

**SKU:** GA541

**Categories:** [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

### ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>Weight Capacity</b>	<a href="#">500 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	7 mph
<b>Operating Range</b>	18 miles
<b>Ground Clearance</b>	4"
<b>Turning Radius</b>	67"
<b>Batteries</b>	NF22
<b>Colors</b>	Red, Silver
<b>Features</b>	Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, PGS 140 Amp Controller
<b>Charger</b>	Off Board

<b>Length</b>	57.5"
<b>Width</b>	24.5"
<b>Height (ground to top of back of seat)</b>	42.5" - 45.5"
<b>Height (ground to top of back of headrest)</b>	50" - 56"
<b>Ground to Top of Seat</b>	25" - 28"
<b>Ground to Top of Deck</b>	6.5"
<b>Number Seat Height Positions</b>	4
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	38.5"
<b>Weight w batteries</b>	265 lbs
<b>Front Half Weight</b>	68 lbs
<b>Rear Half Weight</b>	68 lbs
<b>Seat weight w arms</b>	55 lbs
<b>Batteries Weight</b>	37 lbs ea
<b>Heaviest Piece</b>	68 lbs
<b>Tires</b>	Pneumatic
<b>Front Tire</b>	13"
<b>Rear tire</b>	13"
<b>Casters</b>	N/A
<b>Rear Anti-Tippers</b>	3"
<b>Standard Seat</b>	Captain's Seat
<b>Seat Back Height - no headrest</b>	21" - 24"
<b>Back Height with Headrest</b>	27.5" - 31.5"
<b>Seat Width x Depth</b>	20" X 17"
<b>Seat Color</b>	Gray Vinyl
<b>Warranty Drive Train</b>	2 years



**Warranty Batteries**

1 year