

FDA Class II Medical Device



AVENGER

\$4,699.00

The Golden Avenger—a four-wheel, heavy-duty scooter—provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature’s most challenging terrains.

SKU: GA541

Categories: [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

ADDITIONAL INFORMATION

| | |
|-------------------------|---|
| Medicare Code | Not Coded |
| Weight Capacity | 500 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 7 mph |
| Operating Range | 18 miles |
| Ground Clearance | 4" |
| Turning Radius | 67" |
| Batteries | NF22 |
| Colors | Red, Silver |
| Features | Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, PGS 140 Amp Controller |
| Charger | Off Board |

| | |
|---|----------------|
| Length | 57.5" |
| Width | 24.5" |
| Height (ground to top of back of seat) | 42.5" - 45.5" |
| Height (ground to top of back of headrest) | 50" - 56" |
| Ground to Top of Seat | 25" - 28" |
| Ground to Top of Deck | 6.5" |
| Number Seat Height Positions | 4 |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 38.5" |
| Weight w batteries | 265 lbs |
| Front Half Weight | 68 lbs |
| Rear Half Weight | 68 lbs |
| Seat weight w arms | 55 lbs |
| Batteries Weight | 37 lbs ea |
| Heaviest Piece | 68 lbs |
| Tires | Pneumatic |
| Front Tire | 13" |
| Rear tire | 13" |
| Casters | N/A |
| Rear Anti-Tippers | 3" |
| Standard Seat | Captain's Seat |
| Seat Back Height - no headrest | 21" - 24" |
| Back Height with Headrest | 27.5" - 31.5" |
| Seat Width x Depth | 20" X 17" |
| Seat Color | Gray Vinyl |
| Warranty Drive Train | 2 years |



Warranty Batteries

1 year