



AVENGER

\$4,699.00

The Golden Avenger—a four-wheel, heavy-duty scooter—provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature’s most challenging terrains.

SKU: GA541

Categories: [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range	18 miles
Ground Clearance	4"
Turning Radius	67"
Batteries	NF22
Colors	Red, Silver
Features	Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, PGS 140 Amp Controller
Charger	Off Board

Length	57.5"
Width	24.5"
Height (ground to top of back of seat)	42.5" - 45.5"
Height (ground to top of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
Number Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	38.5"
Weight w batteries	265 lbs
Front Half Weight	68 lbs
Rear Half Weight	68 lbs
Seat weight w arms	55 lbs
Batteries Weight	37 lbs ea
Heaviest Piece	68 lbs
Tires	Pneumatic
Front Tire	13"
Rear tire	13"
Casters	N/A
Rear Anti-Tippers	3"
Standard Seat	Captain's Seat
Seat Back Height - no headrest	21" - 24"
Back Height with Headrest	27.5" - 31.5"
Seat Width x Depth	20" X 17"
Seat Color	Gray Vinyl
Warranty Drive Train	2 years



Warranty Batteries

1 year