

FDA Class II Medical Device



BUZZAROUND CARRYON HD

\$2,649.00

Go on any adventure with this incredibly portable, lightweight travel scooter with a 400 lb. weight capacity. Up to 22 mile coverage with airline-friendly 12.5 AH & spare 6.5 AH lithium ion batteries.

SKU: GB220

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Favorites](#), [Folding](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	4.29 mph
Operating Range (12.5 AH Battery)	14.2 miles
Operating Range (6.5 AH Battery)	7.4 miles
Total Combined Range	21.6 miles
Ground Clearance	4.63"
Ground Clearance to Transaxle/ Motor	2.3"



Turning Radius	48.75"
Type Batteries	12 AH Lithium-Ion, 6.5 AH Lithium-Ion
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	32 Amp
Length	40"
Width	23.5"
Height (ground to top of back of seat)	39"
Maximum Height (folded without seat)	14"
Ground to Top of Seat	23.75"
Ground to Top of Deck	6.6"
Deck to Top of Seat	17"
Front Axle to Rear Axle	28.7"
Weight of Unit (assembled)	77.5 lbs
Main Body Weight	52.5 lbs
Seat with Arms Weight	20 lbs
Batteries Weight	5 lbs each
Heaviest Piece	52.5 lbs
Front Tires	7.25"
Rear Tires	8.25"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	15"
Seat Width x Depth	16" x 18"
Seat Color	Black Vinyl



Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	2 years