

FDA Class II Medical Device



## BUZZAROUND CARRYON HD

**\$2,649.00**

Go on any adventure with this incredibly portable, lightweight travel scooter with a 400 lb. weight capacity. Up to 22 mile coverage with airline-friendly 12.5 AH & spare 6.5 AH lithium ion batteries.

---

**SKU:** GB220

**Categories:** [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Folding](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>User Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4.29 mph
<b>Operating Range (12.5 AH Battery)</b>	14.2 miles
<b>Operating Range (6.5 AH Battery)</b>	7.4 miles
<b>Total Combined Range</b>	21.6 miles
<b>Ground Clearance</b>	4.63"
<b>Ground Clearance to Transaxle/ Motor</b>	2.3"



<b>Turning Radius</b>	48.75"
<b>Type Batteries</b>	12 AH Lithium-Ion, 6.5 AH Lithium-Ion
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	32 Amp
<b>Length</b>	40"
<b>Width</b>	23.5"
<b>Height (ground to top of back of seat)</b>	39"
<b>Maximum Height (folded without seat)</b>	14"
<b>Ground to Top of Seat</b>	23.75"
<b>Ground to Top of Deck</b>	6.6"
<b>Deck to Top of Seat</b>	17"
<b>Front Axle to Rear Axle</b>	28.7"
<b>Weight of Unit (assembled)</b>	77.5 lbs
<b>Main Body Weight</b>	52.5 lbs
<b>Seat with Arms Weight</b>	20 lbs
<b>Batteries Weight</b>	5 lbs each
<b>Heaviest Piece</b>	52.5 lbs
<b>Front Tires</b>	7.25"
<b>Rear Tires</b>	8.25"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Seat Back Height - no headrest</b>	15"
<b>Seat Width x Depth</b>	16" x 18"
<b>Seat Color</b>	Black Vinyl



<b>Warranty Frame</b>	Lifetime
<b>Warranty Drive Train</b>	2 years
<b>Warranty Electronics</b>	2 years
<b>Warranty Batteries</b>	2 years