

BUZZAROUND CARRYON

\$2,811.00

The Buzzaround CarryOn motorized scooter is small in stature and big on fun! This new sporty, super portable scooter will take you where you want to go in comfort and style! It features an infinitely adjustable tiller that can fold flat, a full front and rear lighting package, lithium-ion battery pack (travel battery pack sold separately*), and 300 lb. weight capacity.



SKU: GB120

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Favorites](#), [Folding](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
Weight capacity	300
Drive Wheels	Rear
Maximum speed	3.6 MPH
Operating range	18 MILES
Ground clearance	4"
Ground Clearance to Transaxle/ Motor	1.5"
Turning radius	48.75"

Batteries	15 AH, 6.5 AH (Optional)
Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	38"
Width	19.5"
Height (ground to top of back of seat)	34.5"
Maximum Height (folded without seat)	12"
Ground to Top of Seat	23.5"
Ground to Top of Deck	6"
Deck to Top of Seat	17.5"
Front Axle to Rear Axle	28.7"
Weight w batteries	65.5 LBS
Main Body Weight	47.5 lbs
Seat weight w arms	12 LBS
Front half weight	47.5 LBS
Rear half weight	N/A
Batteries Weight	15 AH 6 lbs, 6.5 AH 3.8 lbs
Heaviest piece	47.5 LBS
Front tire	6"
Rear tire	7"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	13.25"
Seat Width x Depth	15" x 13.75"
Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 Year

Warranty Electronics 2 Years

Warranty Batteries 1 Year