

FDA Class II Medical Device



BUZZAROUND CARRYON

\$2,649.00

Travel anywhere in the Buzzaround CarryOn folding scooter! Whether by plane, cruise or train, this lightweight & portable scooter promises comfort & style on-the-go.

SKU: GB120

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Favorites](#), [Folding](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	3.6 mph
Operating Range (12.5 AH Battery)	14.4 miles
Operating Range (6.5 AH Battery)	7.5 miles
Ground Clearance Under Deck	4"
Ground Clearance to Transaxle/ Motor	1.5"
Turning Radius	48.75"

Type Batteries	12.5 AH Lithium-Ion, 6.5 AH Lithium-Ion
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Length	38"
Width	19.5"
Height (ground to top of back of seat)	34.5"
Ground to Top of Seat	23.5"
Ground to Top of Deck	6"
Deck to Top of Seat	17.5"
Front Axle to Rear Axle	28.7"
Weight of Unit (assembled)	65.5 lbs
Main Body Weight	47.5 lbs
Seat with Arms Weight	12 lbs
Front Half Weight	47.5 lbs
Rear Half Weight	N/A
Batteries Weight	5 lbs each
Heaviest Piece	47.5 lbs
Front Tires	6"
Rear Tires	7"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	13.25"
Seat Width x Depth	15" x 13.75"
Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 years

Warranty Electronics 2 years

Warranty Batteries 2 years