



BUZZAROUND EX 3-WHEEL

\$2,359.00

The Buzzaround EX 3-Wheel features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full-size scooter.

SKU: GB118

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

| | |
|-----------------------------|---|
| Description | EX |
| Medicare Code | Not Coded |
| User Weight Capacity | 350 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 5 mph |
| Operating Range | 18 miles |
| Ground Clearance | 4" |
| Turning Radius | 40.5" |
| Batteries | U1 |
| Features | 90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode |
| Charger | Off Board |
| Length | 45" |



| | |
|--|------------------|
| Width | 22" |
| Height (ground to top of back of seat) | 35"-37" |
| Ground to Top of Seat | 22.75"-24.75" |
| Ground to Top of Deck | 6.75" |
| Deck to Top of Seat | 16"-18" |
| Number Seat Height Positions | 3 |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 34.7" |
| Weight w batteries | 154 lbs |
| Front Half Weight | 46 lbs |
| Rear Half Weight | 36 lbs |
| Seat weight w arms | 23 lbs |
| Batteries Weight | 24 lbs. per pack |
| Heaviest Piece | 46 lbs |
| Tires | Solid |
| Front Tire | 9" |
| Rear tire | 9" |
| Rear Anti-Tippers | 2" |
| Seat Back Height - no headrest | 16" |
| Seat Width x Depth | 18" x 16" |
| Seat Color | Black Vinyl |
| Optional Seat | Stadium Style |
| Optional Seat Back Height - no headrest | 16" |
| Optional Seat Width x Depth | 20" x 18" |
| Optional Seat Color | Black Vinyl |
| Warranty Drive Train | 2 years |



Warranty Batteries

1 year