



BUZZAROUND EX 3-WHEEL

\$2,359.00

The Buzzaround EX 3-Wheel features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full-size scooter.

SKU: GB118

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Description	EX
Medicare Code	Not Coded
User Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	18 miles
Ground Clearance	4"
Turning Radius	40.5"
Type Batteries	2-U1
Standard Colors	Red, Blue
Optional Color Shroud Kits	Grey, Yellow, White
Freewheel Mode	Yes
Electronic Speed Control	Yes



Electro-Mechanical Brakes	Yes
Charger	Off Board
Length	45"
Width	22"
Height (ground to top of back of seat)	35.5"-39.5"
Ground to Top of Seat	22.75"-24.75"
Ground to Top of Deck	6.75"
Deck to Top of Seat	16"-18"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.7"
Weight of Unit (assembled)	154 lbs
Front Half Weight	46 lbs
Rear Half Weight	36 lbs
Seat with Arms Weight	23 lbs
Batteries Weight	24 lbs. per pack
Heaviest Piece	46 lbs
Tires (standard)	Solid, Solid/Low Profile
Front Tires	9"
Rear Tires	9"
Rear Anti-Tippers	2"
Standard Seat	Highback
Seat Back Height - no headrest	18"
Seat Width x Depth	18" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"



Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year