



## BUZZAROUND EX 3-WHEEL

**\$2,540.00**

The Buzzaround EX 3-Wheel Mobility Scooter offers more standard features to make your riding experience better than ever. Featuring a 350 lb. weight capacity, a wraparound Delta Tiller, a new LED battery gauge, a wireless U1 battery, and an ultra-bright LED angle-adjustable headlight, the Buzzaround EX travel scooter is also longer to offer more foot room, with a great turning radius.

---

**SKU:** GB118

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)  
**Tags:** [Best Seller](#)

## ADDITIONAL INFORMATION

<b>Description</b>	EX
<b>Medicare Code</b>	Not Coded
<b>Weight Capacity</b>	<a href="#">350</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 MPH
<b>Operating range</b>	18 MILES
<b>Ground Clearance</b>	4"
<b>Turning radius</b>	40.5"
<b>Batteries</b>	U1



<b>Features</b>	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
<b>Charger</b>	Off Board
<b>Length</b>	45"
<b>Width</b>	22"
<b>Height (ground to top of back of seat)</b>	35"-37"
<b>Ground to Top of Seat</b>	22.75"-24.75"
<b>Ground to Top of Deck</b>	6.75"
<b>Deck to Top of Seat</b>	16"-18"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	34.7"
<b>Weight w batteries</b>	154 LBS
<b>Front Half Weight</b>	46 LBS
<b>Rear Half Weight</b>	36 LBS
<b>Seat weight w arms</b>	23 LBS
<b>Batteries Weight</b>	24 lbs. per pack
<b>Heaviest Piece</b>	46 LBS
<b>Tires</b>	Solid
<b>Front Tire</b>	9"
<b>Rear tire</b>	9"
<b>Rear Anti-Tippers</b>	2"
<b>Seat Back Height - no headrest</b>	16"
<b>Seat Width x Depth</b>	18" x 16"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Stadium Style
<b>Optional Seat Back Height - no headrest</b>	16"



<b>Optional Seat Width x Depth</b>	20" x 18"
<b>Optional Seat Color</b>	Black Vinyl
<b>Warranty Drive Train</b>	2 Year
<b>Warranty Batteries</b>	1 Year