



BUZZAROUND EX 4-WHEEL

\$2,648.00

The Buzzaround EX 4-wheel scooter offers more standard features to make your riding experience better than ever. Featuring a 350 lb. weight capacity, a wrap-around Delta Tiller, a new LED battery gauge, a wireless U1 battery, and a new ultra bright LED angle-adjustable headlight, the new Buzzaround EX mobility scooter is also longer to offer more foot room, with a great turning radius.

SKU: GB148

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Description	EX
Medicare Code	Not Coded
Weight capacity	350
Drive Wheels	Rear
Maximum speed	5 MPH
Operating range	18 MILES
Ground clearance	4"
Turning radius	54"

Batteries	U1
Colors	Blue, Red
Features	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	44"
Width	22"
Height (ground to top of back of seat)	35"-37"
Ground to Top of Seat	22.75"-24.75"
Ground to Top of Deck	6.75"
Deck to Top of Seat	16"-18"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.7"
Weight w batteries	161 LBS
Front half weight	53 LBS
Rear half weight	36 LBS
Seat weight w arms	23 LBS
Batteries Weight	24 lbs. per pack
Heaviest piece	53 LBS
Front tire	9"
Rear tire	9"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	16"
Seat Width x Depth	18" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style

Optional Seat Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 Year
Warranty Electronics	2 Years
Warranty Batteries	1 Year