

BUZZAROUND EX 4-WHEEL

\$2,399.00

The Buzzaround EX 4-Wheel features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full-size scooter.

SKU: GB148

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

Tags: [Best Seller](#)



ADDITIONAL INFORMATION

Description	EX
Medicare Code	Not Coded
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	18 miles
Ground Clearance	4"
Turning Radius	54"
Batteries	U1
Colors	Blue, Red
Features	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board



Length	44"
Width	22"
Height (ground to top of back of seat)	35"-37"
Ground to Top of Seat	22.75"-24.75"
Ground to Top of Deck	6.75"
Deck to Top of Seat	16"-18"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.7"
Weight w batteries	161 lbs
Front Half Weight	53 lbs
Rear Half Weight	36 lbs
Seat weight w arms	23 lbs
Batteries Weight	24 lbs. per pack
Heaviest Piece	53 lbs
Front Tire	9"
Rear tire	9"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	16"
Seat Width x Depth	18" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl

Warranty Drive Train 2 years

Warranty Batteries 1 year