

BUZZAROUND LT 3-WHEEL**\$1,549.00**

The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, this sporty design features a lightweight, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the tiller, a new hard plastic rugged basket, and color through shrouds. Enjoy its tight indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance.

SKU: GB107**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)**Tags:** [New](#)**ADDITIONAL INFORMATION**

Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	7 miles
Ground Clearance	2.5"

Turning Radius	39"
Batteries	2-12AH
Colors	Red
Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	41.75"
Width	19-21"
Height	32.5" - 34"
Ground to Top of Seat	18.75" - 20.25"
Ground to Top of Deck	4.25"
Deck to Top of Seat	14.5" - 16"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	33.5"
Weight w batteries	80 lbs
Front Half Weight	28 lbs
Rear Half Weight	20.5 lbs
Seat weight w arms	13.5 lbs
Batteries Weight	17.5 lbs
Heaviest Piece	28 lbs
Front Tire	8"
Rear tire	8"
Rear Anti-Tippers	2"
Seat Back Height - no headrest	12.75"
Seat Width x Depth	16" x 14"
Seat Color	Black Vinyl
Warranty Drive Train	2 years