

FDA Class II Medical Device



## BUZZAROUND LT 3-WHEEL

**\$1,549.00**

The Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models.

---

**SKU:** GB107

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Description</b>	LT
<b>Medicare Code</b>	K0800
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	7 miles
<b>Ground Clearance Under Deck</b>	2.5"
<b>Turning Radius</b>	39"
<b>Type Batteries</b>	2-12AH
<b>Colors</b>	Red
<b>Freewheel Mode</b>	Yes



<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	50 Amp
<b>Length</b>	41.75"
<b>Width</b>	19-21"
<b>Height (ground to top of back of seat)</b>	31.75"-33.75"
<b>Ground to Top of Seat</b>	19" - 21"
<b>Ground to Top of Deck</b>	4"
<b>Deck to Top of Seat</b>	15"-17"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	.75"
<b>Front Axle to Rear Axle</b>	33.5"
<b>Weight of Unit (assembled)</b>	87 lbs.
<b>Front Half Weight</b>	28 lbs
<b>Rear Half Weight</b>	20.5 lbs
<b>Seat with Arms Weight</b>	13.5 lbs
<b>Batteries Weight</b>	2 packs / 17.5 lbs each
<b>Heaviest Piece</b>	28 lbs
<b>Front Tires</b>	8"
<b>Rear Tires</b>	8"
<b>Rear Anti-Tippers</b>	2"
<b>Seat Back Height - no headrest</b>	12.75"
<b>Seat Width x Depth</b>	16" x 14"
<b>Seat Color</b>	Black Vinyl
<b>Warranty Frame</b>	Lifetime
<b>Warranty Drive Train</b>	2 years



**Warranty Electronics** 2 years

**Warranty Batteries** 1 year