

FDA Class II Medical Device

## **BUZZAROUND LT 3-WHEEL**

**\$1,549.00**



The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, this sporty design features a lightweight, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the tiller, a new hard plastic rugged basket, and color through shrouds. Enjoy its tight indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance.

---

**SKU:** GB107

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

**Tags:** [New](#)

## **ADDITIONAL INFORMATION**

<b>Medicare Code</b>	K0800
<b>Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	7 miles
<b>Ground Clearance</b>	2.5"

<b>Turning Radius</b>	39"
<b>Batteries</b>	2-12AH
<b>Colors</b>	Red
<b>Features</b>	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
<b>Charger</b>	Off Board
<b>Length</b>	41.75"
<b>Width</b>	19-21"
<b>Height</b>	32.5" - 34"
<b>Ground to Top of Seat</b>	18.75" - 20.25"
<b>Ground to Top of Deck</b>	4.25"
<b>Deck to Top of Seat</b>	14.5" - 16"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	.75"
<b>Front Axle to Rear Axle</b>	33.5"
<b>Weight w batteries</b>	80 lbs
<b>Front Half Weight</b>	28 lbs
<b>Rear Half Weight</b>	20.5 lbs
<b>Seat weight w arms</b>	13.5 lbs
<b>Batteries Weight</b>	17.5 lbs
<b>Heaviest Piece</b>	28 lbs
<b>Front Tire</b>	8"
<b>Rear tire</b>	8"
<b>Rear Anti-Tippers</b>	2"
<b>Seat Back Height - no headrest</b>	12.75"
<b>Seat Width x Depth</b>	16" x 14"
<b>Seat Color</b>	Black Vinyl
<b>Warranty Drive Train</b>	2 years