

FDA Class II Medical Device

BUZZAROUND LT 3-WHEEL

\$1,549.00



The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, this sporty design features a lightweight, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the tiller, a new hard plastic rugged basket, and color through shrouds. Enjoy its tight indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance.

SKU: GB107

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	7 miles
Ground Clearance	2.5"



Turning Radius	39"
Type Batteries	2-12AH
Colors	Red
Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	41.75"
Width	19-21"
Height (ground to top of back of seat)	32.5" - 34"
Ground to Top of Seat	18.75" - 20.25"
Ground to Top of Deck	4.25"
Deck to Top of Seat	14.5" - 16"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	33.5"
Weight of Unit (assembled)	80 lbs
Front Half Weight	28 lbs
Rear Half Weight	20.5 lbs
Seat with Arms Weight	13.5 lbs
Batteries Weight	17.5 lbs
Heaviest Piece	28 lbs
Front Tires	8"
Rear Tires	8"
Rear Anti-Tippers	2"
Seat Back Height - no headrest	12.75"
Seat Width x Depth	16" x 14"
Seat Color	Black Vinyl



Warranty Drive Train

2 years