

FDA Class II Medical Device



microTURN

BUZZAROUND LT 4-WHEEL

\$1,699.00

The new Buzzaround LT 4-Wheel model features exclusive **MicroTurn™ technology**, delivering an incredible 41.75" turning radius to navigate even the tightest spaces. It's also the smallest, lightest, and most portable Buzzaround yet—packed with sleek upgrades, a 300 lb. weight capacity, and a smooth indoor/outdoor ride.

SKU: GB108

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	7.4 miles
Ground Clearance Under Deck	2.5"
Turning Radius	41.75"
Type Batteries	2-12AH
Color	Red, Blue



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Length	44"
Width	19-21"
Height (ground to top of seat back)	35" - 36.5"
Ground to Top of Seat	21" - 22.5"
Ground to Top of Deck	4"
Deck to Top of Seat	15"-17"
Number of Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	33.5"
Weight of Unit (assembled)	95 lbs
Front Half Weight	37.5 lbs
Rear Half Weight	20.5 lbs
Seat with Arms Weight	19 lbs
Batteries Weight (pack with 2 batteries)	17.5 lbs
Tires	Solid, Low Profile
Front Tires	8"
Rear Tires	8"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Back Seat Height - no headrest	15"
Seat Width x Depth	17.5" x 15"
Seat Color	Black Vinyl



Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year