

FDA Class II Medical Device



BUZZAROUND LT 4-WHEEL

\$1,699.00

The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, this sporty design features a lightweight, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the tiller, a new hard plastic rugged basket, and color through shrouds. Enjoy its tight indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance.

SKU: GB108

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	7 miles
Ground Clearance	2.5"



Turning Radius	48"
Type Batteries	2-12AH
Color	Red
Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	44"
Width	19-21"
Height (ground to top of seat back)	32.5" - 34"
Ground to Top of Seat	18.75" - 20.25"
Ground to Top of Deck	4.25"
Deck to Top of Seat	14.5" - 16"
Number of Seat Height Positions	3
Size of Increments	.75"
Weight of Unit (assembled)	86 lbs
Front Half Weight	34 lbs
Rear Half Weight	20.5 lbs
Seat with Arms Weight	13.5 lbs
Batteries Weight (pack with 2 batteries)	17.5 lbs
Heaviest Piece	34 lbs
Tires	Solid, Low Profile
Front Tires	8"
Rear Tires	8"
Rear Anti-Tippers	2"
Back Seat Height - no headrest	12.75"
Seat Width x Depth	16" x 14"
Seat Color	Black Vinyl



Warranty Drive Train

2 years