

## BUZZAROUND LX 3-WHEEL

**\$2,648.00**

The Buzzaround LX (Luxury) 3-Wheel Scooter is one of the industry's only "Full-Time" / "All Day" scooters, which is also fully portable for travel.

---

**SKU:** GB119

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)



## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>Weight Capacity</b>	<a href="#">375</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 MPH
<b>Operating range</b>	18 MILES
<b>Ground Clearance</b>	4"
<b>Turning radius</b>	40.5"
<b>Batteries</b>	U1
<b>Colors</b>	Blue, Red
<b>Features</b>	90 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
<b>Charger</b>	Off Board
<b>Length</b>	45"



<b>Width</b>	22"
<b>Height (ground to top of seat back)</b>	37"-39"
<b>Ground to Top of Seat</b>	21.75"-23.75"
<b>Ground to Top of Deck</b>	6.75"
<b>Deck to Top of Seat</b>	15.75"-17.75"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	34.7"
<b>Weight w batteries</b>	168 LBS
<b>Front Half Weight</b>	50 LBS
<b>Rear Half Weight</b>	36 LBS
<b>Seat weight w arms</b>	32 LBS
<b>Batteries Weight</b>	24 lbs./pack X 2= 48lbs
<b>Heaviest Piece</b>	50 LBS
<b>Tires</b>	Solid/Low Profile
<b>Optional Tires</b>	N/A
<b>Front Tire</b>	9"
<b>Rear tire</b>	9"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Highback
<b>Seat Back Height - no headrest</b>	17"
<b>Seat Back Height - with headrest</b>	23"-26"
<b>Seat Width x Depth</b>	17" x 16"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Stadium Style
<b>Optional Seat Back Height - no headrest</b>	16"



<b>Optional Seat Width x Depth</b>	20" x 18"
<b>Optional Seat Color</b>	Black Vinyl
<b>Warranty Drive Train</b>	2 Year
<b>Warranty Batteries</b>	1 Year