

FDA Class II Medical Device



BUZZAROUND LX 3-WHEEL

\$3,099.00

The Buzzaround LX (Luxury) 3-Wheel is a full-size scooter offering comfort and the convenience of true portability to take you where you want to go in style.

SKU: GB119

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	375 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	18 miles
Turning Radius	40.5"
Type Batteries	2-U1
Standard Colors	Red, Blue
Optional Color Shroud Kits	Grey, Yellow, White
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board



Length	45"
Width	22"
Height (ground to top of seat back)	37.5"-39.5"
Ground to Top of Seat	21.75"-23.75"
Ground to Top of Deck	6.75"
Deck to Top of Seat	15.75"-17.75"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.7"
Weight of Unit (assembled)	168 lbs
Front Half Weight	50 lbs
Rear Half Weight	36 lbs
Seat with Arms Weight	32 lbs
Batteries Weight	24 lbs./pack X 2= 48lbs
Heaviest Piece	50 lbs
Tires (standard)	Solid/Low Profile
Optional Tires	N/A
Front Tires	9"
Rear Tires	9"
Rear Anti-Tippers	2"
Standard Seat	Highback
Seat Back Height - no headrest	18"
Seat Back Height - with headrest	23"-26"
Seat Width x Depth	18" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style



Optional Seat Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year