

FDA Class II Medical Device



BUZZAROUND LX 3-WHEEL

\$3,099.00

The Buzzaround LX (Luxury) 3-Wheel is a full-size scooter offering comfort and the convenience of true portability to take you where you want to go in style.

SKU: GB119

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	375 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	18 miles
Ground Clearance	4"
Turning Radius	40.5"
Batteries	U1
Colors	Blue, Red
Features	90 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	45"



Width	22"
Height (ground to top of seat back)	37"-39"
Ground to Top of Seat	21.75"-23.75"
Ground to Top of Deck	6.75"
Deck to Top of Seat	15.75"-17.75"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.7"
Weight w batteries	168 lbs
Front Half Weight	50 lbs
Rear Half Weight	36 lbs
Seat weight w arms	32 lbs
Batteries Weight	24 lbs./pack X 2= 48lbs
Heaviest Piece	50 lbs
Tires	Solid/Low Profile
Optional Tires	N/A
Front Tire	9"
Rear tire	9"
Rear Anti-Tippers	2"
Standard Seat	Highback
Seat Back Height - no headrest	17"
Seat Back Height - with headrest	23"-26"
Seat Width x Depth	17" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"



Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Drive Train	2 years
Warranty Batteries	1 year