



BUZZAROUND XL 3-WHEEL

\$1,799.00

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever. It features a larger wraparound Delta tiller, higher ground clearance of 3.5", and an increased operating range of up to 10 miles on the standard 12 AH batteries.

SKU: GB121-STD

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	10 miles/18 miles
Ground Clearance	3.5"
Turning Radius	37.5"
Type Batteries	2-12AH
Optional Batteries	2-22 AH
Standard Colors	Red, Blue



Optional Color Shroud Kits	White, Yellow, Grey, Tan, Pink, Green
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off-board
Controller	50 Amp
Length	41"
Width	20"
Height (ground to top of back of seat)	33.75" - 35.25"
Ground to Top of Seat	21.75"-23.25"
Deck to Top of Seat	16.5"-18"
Ground to Top of Deck	5"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	31.5"
Weight of Unit (assembled)	93 lbs
Front Half Weight	32 lbs
Rear Half Weight	24 lbs
Seat with Arms Weight	21 lbs
Battery weight	16lbs./30lbs
Heaviest Piece	32 lbs
Tires	Solid, low profile
Front Tires	7.5"
Rear Tires	8"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"



Seat Width x Depth	17"x14.5"
Seat Color	Gray Vinyl
Casters	N/A
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year