

## BUZZAROUND XL 3-WHEEL

**\$1,799.00**

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever. It features a larger wraparound Delta tiller, higher ground clearance of 3.5", and an increased operating range of up to 10 miles on the standard 12 AH batteries.



**SKU:** GB121-STD

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Description</b>	XL
<b>Medicare Code</b>	K0800
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	10 miles/18 miles
<b>Ground Clearance</b>	3.5"
<b>Turning Radius</b>	37.5"
<b>Type Batteries</b>	2-12AH



<b>Optional Batteries</b>	2-22 AH
<b>Standard Colors</b>	Red, Blue
<b>Optional Color Shroud Kits</b>	White, Yellow, Grey, Tan, Pink, Green
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off-board
<b>Controller</b>	50 Amp
<b>Length</b>	41"
<b>Width</b>	20"
<b>Height (ground to top of back of seat)</b>	33.75" - 35.25"
<b>Ground to Top of Seat</b>	21.75"-23.25"
<b>Deck to Top of Seat</b>	16.5"-18"
<b>Ground to Top of Deck</b>	5"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	.75"
<b>Front Axle to Rear Axle</b>	31.5"
<b>Weight of Unit (assembled)</b>	93 lbs
<b>Front Half Weight</b>	32 lbs
<b>Rear Half Weight</b>	24 lbs
<b>Seat with Arms Weight</b>	21 lbs
<b>Battery weight</b>	16lbs./30lbs
<b>Heaviest Piece</b>	32 lbs
<b>Tires</b>	Solid, low profile
<b>Front Tires</b>	7.5"
<b>Rear Tires</b>	8"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style



<b>Seat Back Height - no headrest</b>	14.75"
<b>Seat Width x Depth</b>	17"x14.5"
<b>Seat Color</b>	Gray Vinyl
<b>Casters</b>	N/A
<b>Warranty Frame</b>	Lifetime
<b>Warranty Drive Train</b>	2 years
<b>Warranty Electronics</b>	2 years
<b>Warranty Batteries</b>	1 year