

BUZZAROUND XL 3-WHEEL

\$1,799.00

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever. It features a larger wraparound Delta tiller, higher ground clearance of 3.5", and an increased operating range of up to 10 miles on the standard 12 AH batteries.



SKU: GB121-STD

Categories: [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Description	XL
Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	10 miles
Turning Radius	37.5"
Batteries	12 AH
Colors	Red, Blue

Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off-board
Controller	50 Amp
Length	41"
Width	20"
Floor to Top of Seat	33.75-35.25
Floor to Top of Deck	5
Deck to Top of Seat	16.5"-18"
Number Seat Height Positions	3
Size of Increments	.75"
Size Increments Front Axle to Rear Axle	31.5"
Weight w batteries	93 lbs
Base Weight	N/A
Front Half Weight	32 lbs
Rear Half Weight	24 lbs
Battery weight	16lbs./30lbs
Heaviest Piece	35 lbs
Tires	Solid, low profile
Front Tire	7.5"
Center Tire	N/A
Rear tire	8"
Rear Anti-Tip Wheels	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"
Seat Width x Depth	17"x14.5"



Seat weight w arms	21 lbs
Casters	N/A
Ground Clearance	3.5"
Warranty Drive Train	2 years
Warranty on Electrical Parts	2 years
Warranty Batteries	1 year