



**BUZZAROUND XL 3-WHEEL**

**\$2,107.00**

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever. It features a larger wraparound Delta tiller, higher ground clearance of 3.5", and an increased operating range of up to 10 miles on the standard 12 AH batteries.

---

**SKU:** GB121B-STD

**Categories:** [3-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

**Tags:** [New](#)

**ADDITIONAL INFORMATION**

<b>Description</b>	XL
<b>Medicare Code</b>	K0800
<b>Weight capacity</b>	<a href="#">300</a>
<b>Drive Wheels</b>	Rear
<b>Maximum speed</b>	4 MPH
<b>Operating range</b>	10 MILES
<b>Turning radius</b>	37.5"
<b>Batteries</b>	12 AH
<b>Colors</b>	Red, Blue

<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off-board
<b>Controller</b>	50 Amp
<b>Length</b>	41"
<b>Width</b>	20"
<b>Floor to top of seat</b>	33.75-35.25
<b>Floor to Top of Deck</b>	5
<b>Deck to Top of Seat</b>	16.5"-18"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	.75"
<b>Size Increments Front Axle to Rear Axle</b>	31.5"
<b>Weight w batteries</b>	93 LBS
<b>Base weight</b>	N/A
<b>Front half weight</b>	32 LBS
<b>Rear half weight</b>	24 LBS
<b>Battery weight</b>	16lbs./30lbs
<b>Heaviest piece</b>	35 LBS
<b>Tires</b>	Solid, low profile
<b>Front tire</b>	7.5"
<b>Center tire</b>	N/A
<b>Rear tire</b>	8"
<b>Rear Anti-Tip Wheels</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Seat Back Height - no headrest</b>	14.75"
<b>Seat Width x Depth</b>	17"x14.5"



<b>Seat weight w arms</b>	21 LBS
<b>Casters</b>	N/A
<b>Ground clearance</b>	3.5"
<b>Warranty Frame</b>	Lifetime
<b>Warranty Drive Train</b>	2 Year
<b>Warranty electrical parts</b>	2 Year
<b>Warranty Batteries</b>	1 Year