

BUZZAROUND XL 4-WHEEL

\$1,999.00

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever.

SKU: GB124-STD

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)



ADDITIONAL INFORMATION

Description	XL
Medicare Code	K0800
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	10 miles/18 miles
Ground Clearance Under Deck	3.5"
Turning Radius	45.5"
Type Batteries	2-12AH
Optional Batteries	2-22 AH
Colors	Blue, Red



Optional Color Shroud Kits	White, Yellow, Grey, Tan, Pink, Green
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Length	40"
Width	20"
Height (ground to top of back of seat)	33.75"-35.25"
Ground to Top of Seat	21.75"-23.25"
Ground to Top of Deck	5"
Deck to Top of Seat	16.5"-18"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	30.5"
Weight of Unit (assembled)	98 lbs
Front Half Weight	37 lbs
Rear Half Weight	24 lbs
Seat with Arms Weight	21 lbs
Heaviest Piece	37 lbs
Front Tires	7.5"
Rear Tires	8"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"
Seat Width x Depth	17" x 14.5"
Seat Color	Black Vinyl



Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year