

## BUZZAROUND XL 4-WHEEL

**\$1,999.00**

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever.

---

**SKU:** GB124-STD

**Categories:** [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

**Tags:** [New](#)



### ADDITIONAL INFORMATION

<b>Description</b>	XL
<b>Medicare Code</b>	K0800
<b>Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	10 miles
<b>Ground Clearance</b>	3.5"
<b>Turning Radius</b>	45.5"
<b>Batteries</b>	12 AH
<b>Colors</b>	Blue, Red
<b>Features</b>	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode

<b>Charger</b>	Off Board
<b>Length</b>	40"
<b>Width</b>	20"
<b>Height (ground to top of back of seat)</b>	33.75"-35.25"
<b>Ground to Top of Seat</b>	21.75"-23.25"
<b>Ground to Top of Deck</b>	5"
<b>Deck to Top of Seat</b>	16.5"-18"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	.75"
<b>Front Axle to Rear Axle</b>	31.5"
<b>Weight w batteries</b>	98 lbs
<b>Front Half Weight</b>	37 lbs
<b>Rear Half Weight</b>	24 lbs
<b>Seat weight w arms</b>	21 lbs
<b>Heaviest Piece</b>	37 lbs
<b>Front Tire</b>	7.5"
<b>Rear tire</b>	8"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Seat Back Height - no headrest</b>	14.75"
<b>Seat Width x Depth</b>	17" x 14.5"
<b>Seat Color</b>	Black Vinyl
<b>Warranty Drive Train</b>	2 years
<b>Warranty Batteries</b>	1 year