

BUZZAROUND XL 4-WHEEL

\$1,999.00

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever.

SKU: GB124-STD

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)



ADDITIONAL INFORMATION

Description	XL
Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	10 miles
Ground Clearance	3.5"
Turning Radius	45.5"
Batteries	12 AH
Colors	Blue, Red
Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode

Charger	Off Board
Length	40"
Width	20"
Height (ground to top of back of seat)	33.75"-35.25"
Ground to Top of Seat	21.75"-23.25"
Ground to Top of Deck	5"
Deck to Top of Seat	16.5"-18"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	31.5"
Weight w batteries	98 lbs
Front Half Weight	37 lbs
Rear Half Weight	24 lbs
Seat weight w arms	21 lbs
Heaviest Piece	37 lbs
Front Tire	7.5"
Rear tire	8"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"
Seat Width x Depth	17" x 14.5"
Seat Color	Black Vinyl
Warranty Drive Train	2 years
Warranty Batteries	1 year