



BUZZAROUND XLS-HD 4-WHEEL

\$2,099.00

This popular Buzzaround XLSHD scooter features compact portability with a heavy-duty weight capacity. It also features our Comfort Spring Suspension for a smoother ride. Standard colors are Red and Blue and optional color shrouds include Green, Tan, White, Yellow, Pink and Grey.

SKU: GB124-SHZ

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Description	XLSHD
Medicare Code	K0801
Weight Capacity	325 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	9 miles
Ground Clearance Under Deck	4"
Turning Radius	44"
Batteries	12 AH
Colors	Blue, Red

Features	50 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Charger	Off Board
Length	40"
Width	20"
Height (ground to top of back of seat)	34.25"-35.75"
Ground to Top of Seat	22.25"-23.75"
Ground to Top of Deck	5.5"
Deck to Top of Seat	16.5"-18"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	31.5"
Weight w batteries	105 lbs
Front Half Weight	40 lbs
Rear Half Weight	28 lbs
Seat weight w arms	21 lbs
Batteries Weight	16 lbs. / 30 lbs
Heaviest Piece	40 lbs
Rear Anti-Tippers	2"
Front Tire	7.5"
Rear tire	8"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"
Seat Width x Depth	17" x 14.5"
Seat Color	Black Vinyl
Warranty Drive Train	2 years
Warranty Batteries	1 year