



BUZZAROUND XLS-HD 4-WHEEL

\$2,099.00

This popular Buzzaround XLSHD scooter features compact portability with a heavy-duty weight capacity. It also features our Comfort Spring Suspension for a smoother ride. Standard colors are Red and Blue and optional color shrouds include Green, Tan, White, Yellow, Pink and Grey.

SKU: GB124-SHZ

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Description	XLSHD
Medicare Code	K0801
User Weight Capacity	325 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	9 miles/18 miles
Ground Clearance Under Deck	4"
Turning Radius	44"
Type Batteries	2-12AH
Optional Batteries	2-22 AH



Colors	Blue, Red
Optional Color Shroud Kits	White, Yellow, Grey, Tan, Pink, Green
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Length	40"
Width	20"
Height (ground to top of back of seat)	34.25"-35.75"
Ground to Top of Seat	22.25"-23.75"
Ground to Top of Deck	5.5"
Deck to Top of Seat	16.5"-18"
Number Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	30.5"
Weight of Unit (assembled)	105 lbs
Front Half Weight	40 lbs
Rear Half Weight	28 lbs
Seat with Arms Weight	21 lbs
Batteries Weight	16 lbs. / 30 lbs
Heaviest Piece	40 lbs
Front Tires	7.5"
Rear Tires	8"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	14.75"



Seat Width x Depth	17" x 14.5"
Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year