



BUZZAROUND XLS-HD 4-WHEEL

\$2,099.00

This popular Buzzaround XLSHD scooter features compact portability with a heavy-duty weight capacity. It also features our Comfort Spring Suspension for a smoother ride. Standard colors are Red and Blue and optional color shrouds include Green, Tan, White, Yellow, Pink and Grey.

SKU: GB124-SHZ

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

| | |
|------------------------------------|-------------------------|
| Description | XLSHD |
| Medicare Code | K0801 |
| User Weight Capacity | 325 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 4 mph |
| Operating Range | 9 miles/18 miles |
| Ground Clearance Under Deck | 4" |
| Turning Radius | 44" |
| Type Batteries | 2-12AH |
| Optional Batteries | 2-22 AH |



| | |
|-----------------------------------------------|---------------------------------------|
| Colors | Blue, Red |
| Optional Color Shroud Kits | White, Yellow, Grey, Tan, Pink, Green |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Charger | Off Board |
| Controller | 50 Amp |
| Length | 40" |
| Width | 20" |
| Height (ground to top of back of seat) | 34.25"-35.75" |
| Ground to Top of Seat | 22.25"-23.75" |
| Ground to Top of Deck | 5.5" |
| Deck to Top of Seat | 16.5"-18" |
| Number Seat Height Positions | 3 |
| Size of Increments | .75" |
| Front Axle to Rear Axle | 30.5" |
| Weight of Unit (assembled) | 105 lbs |
| Front Half Weight | 40 lbs |
| Rear Half Weight | 28 lbs |
| Seat with Arms Weight | 21 lbs |
| Batteries Weight | 16 lbs. / 30 lbs |
| Heaviest Piece | 40 lbs |
| Front Tires | 7.5" |
| Rear Tires | 8" |
| Rear Anti-Tippers | 2" |
| Standard Seat | Stadium Style |
| Seat Back Height - no headrest | 14.75" |



| | |
|-----------------------------|-------------|
| Seat Width x Depth | 17" x 14.5" |
| Seat Color | Black Vinyl |
| Warranty Frame | Lifetime |
| Warranty Drive Train | 2 years |
| Warranty Electronics | 2 years |
| Warranty Batteries | 1 year |