

FDA Class II Medical Device



BUZZAROUND XL+ 4-WHEEL

\$2,349.00

The Buzzaround XL+ 4-Wheel Mobility Scooter offers compact portability with a 325 lb. weight capacity, smooth comfort suspension, and MicroTurn™ Technology for an ultra-tight 39" turning radius. With a lightweight lithium battery, up to a 9.2-mile range, and modern comfort features, it delivers reliable performance and easy maneuverability wherever you go.

SKU: GB152-XL+

Categories: [4-Wheel Scooters](#), [Buzzaround Portable Mobility Scooters](#), [Mobility Scooters](#), [Portable](#), [Travel Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Description	XL+
Medicare Code	K0801
User Weight Capacity	325 lbs
Drive Wheels	Rear
Maximum Speed	4.7 mph
Operating Range	9.2/18 miles
Ground Clearance Under Deck	3.5"
Turning Radius	39"

Type Batteries	12AH 24V Lithium-Ion
Colors	Blue, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Length	42"
Width	22"
Height (ground to top of back of seat)	35.5"-37.5"
Ground to Top of Seat	21.5" - 23.5"
Ground to Top of Deck	6"
Deck to Top of Seat	14.5"-16.5"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30.6"
Weight of Unit (assembled)	97 lbs.
Front Half Weight	42.5 lbs.
Rear Half Weight	28 lbs
Seat with Arms Weight	16 lbs.
Batteries Weight	7 lbs.
Front Tires	8"
Rear Tires	8"
Rear Anti-Tippers	2"
Seat Back Height - no headrest	17"
Seat Width x Depth	16"W-15.5"D
Seat Color	Black

Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year