



COMPANION (3-WHEEL) FULL SIZE

\$3,149.00

Golden's popular Companion 3-wheel scooter has a sporty new look and upgraded comfort features! This next-generation model is the ultimate value in a full-time, yet portable scooter.

SKU: GC340

Categories: [3-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

ADDITIONAL INFORMATION

Medicare Code	K0807
User Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	19 miles
Ground Clearance	4.5"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	45"
Type Batteries	2-U1/ up to 50AH
Colors	Blue, Grey, Red
Charger	Portable 3.5 Amp Super Charger Off Board



Length	46"
Width	24.5"
Ground to Top of Seat	23" - 25"
Ground to Top of Deck	6.5"
Deck to Top of Seat	16.5" - 18.5"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.5"
Weight of Unit (assembled)	207 lbs
Tires (standard)	Flat Free / Non Marking Solid
Front Tires	9"
Rear Tires	10.75"
Seat Back Height - no headrest	18"
Seat Width x Depth	18" x 16"
Seat Color	Black
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black
Seat Options	Pneumatic Rise Seat, Power Elevating Seat - 6" with Lockout
Warranty Drive Train	2 years
Warranty Batteries	1 year