

FDA Class II Medical Device



COMPANION (3-WHEEL) MID SIZE

\$2,759.00

The Companion GC240 scooter will take you the places you've had trouble getting to. Ideal for trips to the mall, dining out or an afternoon at the park, this midsize scooter provides great style with high performance and easy disassembly.

SKU: GC240

Categories: [3-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

ADDITIONAL INFORMATION

Medicare Code	K0807
User Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	19 miles
Ground Clearance	3.75"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	37"
Type Batteries	U1
Colors	Blue, Red



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Portable 3.5 Amp Off Board
Controller	90 Amp
Length	42.25"
Width	24"
Height (ground to top of back of seat)	37.5" - 39.5"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	22.25" - 24.30"
Ground to Top of Deck	6.25"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight of Unit (assembled)	170.8 lbs
Batteries Weight	21 lbs ea
Tires (standard)	Solid/Low Profile
Front Tires	9"
Rear Tires	10.75"
Rear Anti-Tippers	3"
Standard Seat	Quilted High Back
Seat Back Height - no headrest	16"
Seat Width x Depth	18" x 16"
Seat Color	Black
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	18"



Optional Seat Width x Depth	18"x 16"
Optional Seat Color	Black
Seat Options	Pneumatic Rise Seat
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year