

FDA Class II Medical Device

COMPANION (3-WHEEL) MID SIZE

\$2,759.00

The Companion GC240 scooter will take you the places you've had trouble getting to. Ideal for trips to the mall, dining out or an afternoon at the park, this midsize scooter provides great style with high performance and easy disassembly.



SKU: GC240

Categories: [3-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

ADDITIONAL INFORMATION

Medicare Code	K0807
User Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	15 miles
Ground Clearance	3.75"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	37"
Batteries	U1
Colors	Blue, Red



Charger	Portable 3.5 Amp Off Board
Length	42.25"
Width	24"
Ground to Top of Seat	22.25" - 24.30"
Ground to Top of Deck	6.25"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight w batteries	170.8 lbs
Batteries Weight	21 lbs ea
Heaviest Piece	48 lbs
Tires	Flat Free / Non Marking Solid
Front Tire	9"
Rear tire	10.75"
Rear Anti-Tippers	3"
Standard Seat	Quilted High Back
Seat Back Height - no headrest	18"
Seat Width x Depth	18" x 16"
Seat Color	Black
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"
Optional Seat Width x Depth	20"x 18"
Optional Seat Color	Black
Seat Options	Pneumatic Rise Seat
Warranty Drive Train	2 years
Warranty Batteries	1 year