

COMPANION (4-WHEEL) FULL SIZE

\$3,299.00

Golden's popular Companion 4-wheel scooter has a sporty new look and upgraded comfort features! This next-generation model is rugged, dependable and an outstanding value for a full-time, yet portable 4-wheel scooter.



SKU: GC440

Categories: [4-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

ADDITIONAL INFORMATION

Medicare Code	K0807
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	19 miles
Ground Clearance	4.25"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	56"
Batteries	2-U1/ up to 50AH
Colors	Blue, Grey, Red

Charger	Portable 3.5 Amp Super Charger Off Board
Length	48"
Width	24.5"
Ground to Top of Seat	23.5" - 25.5"
Ground to Top of Deck	6.5"
Deck to Top of Seat	17" - 19"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	35.5"
Weight w batteries	222 lbs
Tires	Flat Free / Non Marking Solid
Front Tire	10.75"
Rear tire	10.75"
Rear Anti-Tippers	2"
Seat Back Height - no headrest	18"
Seat Width x Depth	18" x 16"
Seat Color	Black
Optional Seat	Stadium Style
Optional Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black
Seat Options	Pneumatic Rise Seat, Power Elevating Seat - 6" with Lockout
Warranty Drive Train	2 years
Warranty Batteries	1 year