



## COMPANION (4-WHEEL) FULL SIZE

**\$3,299.00**

Golden's popular Companion 4-wheel scooter has a sporty new look and upgraded comfort features! This next-generation model is rugged, dependable and an outstanding value for a full-time, yet portable 4-wheel scooter.

---

**SKU:** GC440

**Categories:** [4-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	K0807
<b>Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 mph
<b>Operating Range</b>	19 miles
<b>Ground Clearance</b>	4.25"
<b>Ground Clearance to Transaxle/Motor</b>	2.5"
<b>Turning Radius</b>	56"
<b>Batteries</b>	2-U1/ up to 50AH
<b>Colors</b>	Blue, Grey, Red



<b>Charger</b>	Portable 3.5 Amp Super Charger Off Board
<b>Length</b>	48"
<b>Width</b>	24.5"
<b>Ground to Top of Seat</b>	23.5" - 25.5"
<b>Ground to Top of Deck</b>	6.5"
<b>Deck to Top of Seat</b>	17" - 19"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	35.5"
<b>Weight w batteries</b>	222 lbs
<b>Tires</b>	Flat Free / Non Marking Solid
<b>Front Tire</b>	10.75"
<b>Rear tire</b>	10.75"
<b>Rear Anti-Tippers</b>	2"
<b>Seat Back Height - no headrest</b>	18"
<b>Seat Width x Depth</b>	18" x 16"
<b>Seat Color</b>	Black
<b>Optional Seat</b>	Stadium Style
<b>Optional Back Height - no headrest</b>	16"
<b>Optional Seat Width x Depth</b>	20" x 18"
<b>Optional Seat Color</b>	Black
<b>Seat Options</b>	Pneumatic Rise Seat, Power Elevating Seat - 6" with Lockout
<b>Warranty Drive Train</b>	2 years
<b>Warranty Batteries</b>	1 year