

FDA Class II Medical Device



COMPANION (4-WHEEL) FULL SIZE

\$3,299.00

Golden's popular Companion 4-wheel scooter has a sporty new look and upgraded comfort features! This next-generation model is rugged, dependable and an outstanding value for a full-time, yet portable 4-wheel scooter.

SKU: GC440

Categories: [4-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Mobility Scooters](#), [Portable](#)

ADDITIONAL INFORMATION

Medicare Code	K0807
User Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	19 miles
Ground Clearance	4.25"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	54"
Type Batteries	2-U1/ up to 50AH
Colors	Blue, Grey, Red



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Portable 3.5 Amp Super Charger Off Board
Controller	90 Amp
Length	48"
Width	24"
Height (ground to top of back of seat)	38" - 40"
Ground to Top of Seat	23" - 25"
Ground to Top of Deck	6.5"
Deck to Top of Seat	17" - 19"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	35.5"
Weight of Unit (assembled)	222 lbs
Tires (standard)	Solid/Low Profile
Front Tires	10.75"
Rear Tires	10.75"
Rear Anti-Tippers	2"
Standard Seat	Highback
Seat Back Height - no headrest	18"
Seat Width x Depth	18" x 16"
Seat Color	Black
Optional Seat	Stadium Style
Optional Back Height - no headrest	16"
Optional Seat Width x Depth	20" x 18"



Optional Seat Color	Black
Seat Options	Pneumatic Rise Seat, Power Elevating Seat – 6" with Lockout
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year