

FDA Class II Medical Device



COMPANION HD

\$3,799.00

The new Golden Companion Heavy-Duty 3-wheel scooter combines a modern design with high performance features to offer the ultimate indoor or outdoor driving experience at a weight capacity of 500 lbs.

FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with physical disabilities.

SKU: GC540

Categories: [3-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	K0808
User Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	6 mph
Operating Range	16 miles
Ground Clearance to Center Deck	6.25"
Ground Clearance to Transaxle/Motor	3.75"



Turning Radius	51"
Type Batteries	2-NF22
Optional Batteries	N/A
Colors	Blue, Grey, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Portable 8 Amp Off Board
Controller	120 Amp
Length	48"
Width	30.5"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	24.5" - 27.5"
Ground to Top of Deck	8"
Deck to Top of Seat	16.5" - 19.5"
Number Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	36"
Weight of Unit (assembled)	225.5 lbs
Tires (standard)	Flat Free / Non Marking Solid
Front Tires	10.75"
Rear Tires	10.75"
Front & Rear Anti-Tippers	3"
Standard Seat	Extra Wide High Back
Seat Back Height - no headrest	18"
Seat Width x Depth	22" x 20"
Seat Color	Black



Seat Rotation	90 Degree either direction 360 Degrees
Optional Seat	Captain's Chair
Optional Seat Back Height - no headrest	18"
Optional Seat Width x Depth	20" x 20"
Optional Seat Color	Black
Power Elevating Seat - 6" with Lockout	Yes
Warranty on Frame	5 Year
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year