



COMPANION HD

\$4,799.00

The new Golden Companion Heavy-Duty 3-wheel scooter combines a modern design with high performance features to offer the ultimate indoor or outdoor driving experience at a weight capacity of 500 lbs.

SKU: GC540

Categories: [3-Wheel Scooters](#), [Companion Full-Size Mobility Scooters](#), [Favorites](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	K0808
Weight Capacity	500
Drive Wheels	Rear
Maximum Speed	6 MPH
Operating range	16 MILES
Ground Clearance to Center Deck	6.25"
Ground Clearance to Transaxle/Motor	3.75"
Turning radius	51"
Batteries	2-NF22
Air Travel Safe	Yes



Colors	Blue, Grey, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Energy Efficient Auto Turn Off (Sleep Mode)	Yes
Charger	Portable 8 Amp Off Board
Controller	120 Amp
Length	48"
Width	26"
Ground to Top of Seat	24.5" - 27.5"
Ground to Top of Deck	8"
Deck to Top of Seat	16.5" - 19.5"
Number Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	36"
Weight of Unit (assembled)	225.5 lbs
Tires	Flat Free / Non Marking Solid
Front Tire	10.25"
Rear tire	10.75"
Front & Rear Anti-Tippers	Yes
Standard Seat	Extra Wide High Back
Seat Back Height - no headrest	18"
Seat Width x Depth	22" x 20"
Seat Color	Black
Transportation Seat Lock	Yes, Self Locking
Seat Rotation	90 Degree either direction 360 Degrees
Optional Seat	Captain's Chair



Optional Seat Back Height - no headrest	18"
Optional Seat Color	Black
Power Elevating Seat - 6" with Lockout	Yes
Warranty on Frame	5 Year
Warranty Drive Train	2 Year
Warranty on Electrical Parts	2 YR
Warranty Batteries	1 Year