



EAGLE

\$4,599.00

If you're looking for the ultimate outdoor driving experience, look no further than the new Golden Eagle™.

SKU: GR595

Categories: [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range	23 miles
Ground Clearance	5"
Turning Radius	69"
Batteries	2-75 AH, NF22
Colors	Red
Features	140 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, Slide Seat Bracket
Charger	Off Board
Length	56.5"

Width	27"
Height (ground to top of back of seat)	49" - 50"
Height (ground to top of back of headrest)	57" - 62.5"
Ground to Top of Seat	27" - 28"
Ground to Top of Deck	7.25"
Deck to Top of Seat	19.75" - 20.75"
Number Seat Height Positions	2
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight w batteries	302.5 lbs
Base Weight	152 lbs
Seat weight w arms	49.5 lbs
Batteries Weight	37 lbs ea (NF22), 52 lbs ea (75 AH)
Heaviest Piece	179 lbs
Tires	Pneumatic
Front Tire	13"
Rear tire	13"
Standard Seat	Captain's High-Back
Seat Back Height - no headrest	23.75"
Seat Back Height with headrest	30.75" - 36.25"
Seat Width x Depth	20" x 18"
Seat Color	Black
Warranty Drive Train	2 years
Warranty Batteries	1 year