



EAGLE

\$4,599.00

If you're looking for the ultimate outdoor driving experience, look no further than the new Golden Eagle™.

FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with physical disabilities.

SKU: GR595

Categories: [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range	23 miles
Ground Clearance	5"
Ground Clearance to Transaxle/Motor	3.5"
Turning Radius	69"
Type Batteries	2-75 AH, NF22
Colors	Grey, Red, Tan



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
Rear Anti-Tippers	3"
Length	56.5"
Width	27"
Height (ground to top of back of seat)	49" - 50"
Height (ground to top of back of headrest)	57" - 62.5"
Ground to Top of Seat	27" - 28"
Ground to Top of Deck	7.25"
Deck to Top of Seat	19.75" - 20.75"
Number Seat Height Positions	2
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit (assembled)	302.5 lbs
Base Weight	179 lbs.
Seat with Arms Weight	49.5 lbs
Batteries Weight	37 lbs ea (NF22), 52 lbs ea (75 AH)
Tires (standard)	Pneumatic
Front Tires	13"
Rear Tires	13"
Standard Seat	Captain's High-Back
Seat Back Height - no headrest	23.75"
Seat Back Height with headrest	30.75" - 36.25"



Seat Width x Depth	20" x 18"
Seat Color	Black
Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year