



## EAGLE

**\$4,599.00**

If you're looking for the ultimate outdoor driving experience, look no further than the new Golden Eagle™.

---

**SKU:** GR595

**Categories:** [4-Wheel Scooters](#), [Full-Size Scooters](#), [Heavy Duty Scooters](#), [Mobility Scooters](#), [Off-Road Mobility Scooters](#)

**Tags:** [New](#)

### ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>Weight Capacity</b>	<a href="#">500 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	8 mph
<b>Operating Range</b>	23 miles
<b>Ground Clearance</b>	5"
<b>Turning Radius</b>	69"
<b>Batteries</b>	2-75 AH, NF22
<b>Colors</b>	Red
<b>Features</b>	140 Amp Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode, Slide Seat Bracket
<b>Charger</b>	Off Board
<b>Length</b>	56.5"

<b>Width</b>	27"
<b>Height (ground to top of back of seat)</b>	49" - 50"
<b>Height (ground to top of back of headrest)</b>	57" - 62.5"
<b>Ground to Top of Seat</b>	27" - 28"
<b>Ground to Top of Deck</b>	7.25"
<b>Deck to Top of Seat</b>	19.75" - 20.75"
<b>Number Seat Height Positions</b>	2
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	40"
<b>Weight w batteries</b>	302.5 lbs
<b>Base Weight</b>	152 lbs
<b>Seat weight w arms</b>	49.5 lbs
<b>Batteries Weight</b>	37 lbs ea (NF22), 52 lbs ea (75 AH)
<b>Heaviest Piece</b>	179 lbs
<b>Tires</b>	Pneumatic
<b>Front Tire</b>	13"
<b>Rear tire</b>	13"
<b>Standard Seat</b>	Captain's High-Back
<b>Seat Back Height - no headrest</b>	23.75"
<b>Seat Back Height with headrest</b>	30.75" - 36.25"
<b>Seat Width x Depth</b>	20" x 18"
<b>Seat Color</b>	Black
<b>Warranty Drive Train</b>	2 years
<b>Warranty Batteries</b>	1 year