

FDA Class II Medical Device



LITERIDER 3-WHEEL

\$2,099.00

The LiteRider™ three-wheel scooter features up to a 12-mile range, a Delta Tiller standard, more attractive styling, and more vibrant colors. A modern tubular frame design supports its 300 lb. weight capacity and its completely wireless disassembly makes transporting a breeze.

SKU: GL111

Categories: [3-Wheel Scooters](#), [LiteRider](#)
[Portable Mobility Scooters](#), [Mobility Scooters](#),
[Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	12 miles
Ground Clearance Under Deck	3.75"
Ground Clearance to Transaxle/ Motor	2.5"
Turning Radius	39"
Type Batteries	2-22 AH



Colors	Blue, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp
Length	41.5"
Width	22"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	23.5" - 25.5"
Ground to Top of Deck	6"
Deck to Top of Seat	18" - 20"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight of Unit (assembled)	125 lbs
Front Half Weight	44 lbs
Rear Half Weight	28 lbs
Seat with Arms Weight	21 lbs
Batteries Weight	29 lbs/pack
Heaviest Piece	44 lbs
Front Tires	9"
Rear Tires	9"
Rear Anti-Tippers	2"
Standard Seat	Highback
Seat Back Height - no headrest	16"
Seat Width x Depth	18" x 16"



Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Frame	Lifetime
Warranty Drive Train	2 years
Warranty Electronics	2 years
Warranty Batteries	1 year